

INTELLEGO

TECHNOLOGIES

TECHNOLOGY SAVING LIFE AND MONEY

Dosimeter technology background

- Been around since 1800s
- Different types:
 - Radiometer
 - Chemical indicators
 - Biological indicators
 - Etc.
- Used in research, industry, healthcare,
- etc.



Why use a chemical colour changing dosimeter with UV in healthcare?

- Colour changing dosimeters already being used in healthcare:
 - Autoclaves
 - Steam
 - Heat
 - Etc.



US market research (110 hospitals):

> 73% says they would use a dosimeter (some hospitals did not use UV and/or does not plan to use it)

Machine specific for largest actors in industry – on average 80% said they would use the indicators

With 2 of the top 5 companies – 100% of the hospitals said they would use a dosimeter

Colour changing colour indicator most popular way to continuously measure UVC for every cycle

What part does a colour changing dosimeter have in a standard for UVC disinfection ?

- Reference to other standards:
 - Colour changing dosimeter already used in healthcare
 - Recommended by CDC and FDA
- Easy way to monitor UVC exposure for every UVC cycle
- Colour changing dosimeter wont replace radiometer, microbiology testing



CDC Centers for Disease Control and Prevention
CDC 24/7: Saving Lives, Protecting People™

Search

Infection Control

Infection Control > Disinfection and Sterilization

- Disinfection and Sterilization
- Updates
- Authors
- Executive Summary
- Introduction, Methods, Definition of Terms
- A Rational Approach to Disinfection and Sterilization
- Disinfection of Healthcare Equipment
- Factors Affecting the Efficacy of Disinfection and Sterilization
- Cleaning

Performance Indicators

Guideline for Disinfection and Sterilization in Healthcare Facilities (2008)

1. Monitor adherence to high-level disinfection and/or sterilization guidelines for endoscopes on a regular basis. This monitoring should include ensuring the proper training of persons performing reprocessing and their adherence to all endoscope reprocessing steps, as demonstrated by competency testing at commencement of employment and annually.
2. Develop a mechanism for the occupational health service to report all adverse health events potentially resulting from exposure to disinfectants and sterilants; review such exposures; and implement engineering, work practice, and PPE to prevent future exposures.
3. Monitor possible sterilization failures that resulted in instrument recall. Assess whether additional training of personnel or equipment maintenance is required.

[Top of Page](#)

Page last reviewed: September 18, 2016

FDA U.S. FOOD & DRUG ADMINISTRATION

Home / Sterilization Process Controls

Sterilization Process Controls

Share Tweet LinkedIn Email Print

Production and Process Controls Subsystem Sterilization Process Controls

Parent Section N/A

- [Inspectional Objectives](#)
- [Decision Flow Chart](#)

Colour changing chemical indicators development - UVC

- Start with what the market wants (Ips, EVS, etc):
 - Reliable
 - Easy to use
- Photochromic ink
 - Chemical reaction leading to colour change
- Major development points:
 - Wavelength specificity – visible light
 - Stability
 - Repeatability
 - Angles – vertical vs horizontal
 - Temperature, humidity
 - Safety

How can a colour changing indicators help?

- Increase understanding for UVC
- Help to identify shadowed areas
 - 30-50% of surfaces in shadowed areas during cycle*
- Areas further away
- Decrease disinfection time with up to 75%
- **Potential to increase disinfection levels**

*John M. Boyce, MD;1,2 Patricia A. Farrel, MT et.al, “Impact of Room Location on UV-C Irradiance and UV-C Dosage and Antimicrobial Effect Delivered by a Mobile UV-C Light Device”, infection control & hospital epidemiology june 2016, vol. 37, no. 6

Marie Lindblad, Fredrik Huss, MD et.al, “Ultraviolet-C decontamination of a hospital room: Amount of UV light needed”, JBUR 5947 No. of Pages

Does it work?


Ultraviolet-C (UV-C) monitoring *made ridiculously simple:* UV-C dose indicators for convenient measurement of UV-C dosing

Jennifer L. Cadnum, BS¹; Annette L. Jencson, CIC²; Sarah Redmond, MD³; Thriveen Sankar Chittoor Mana, MS³
& Curtis J. Donskey, MD¹⁻³
1. Research Service, Louis Stokes Cleveland Veterans Affairs Medical Center, Cleveland OH,
2. The Cleveland VA Medical Research and Education Foundation, Cleveland, OH
3. Department of Medicine, Division of Infectious Diseases, Case Western Reserve University, Cleveland, OH

Poster #1215

Background

- Ultraviolet-C (UV-C) light is increasingly used as an adjunct to standard cleaning in healthcare facilities
- Most facilities do not have a means to measure UV-C to determine if effective doses are being delivered
- We tested the efficacy of 2 easy-to-use colorimetric indicators for monitoring UV-C dosing in comparison to log reductions in pathogens




Methods

- In a laboratory setting, we exposed methicillin-resistant *Staphylococcus aureus* (MRSA) and *Clostridioides difficile* spores on steel disk carriers to UV-C for varying fluence exposures ranging from 10,000 to 100,000 µJ/cm²
- The UV-C indicators were placed adjacent to the carriers
- Change in color of the indicators was correlated with dose and log₁₀ CFU reductions

Results

Figure 1: Reduction of *C. difficile* (A) and MRSA (B) with varying fluence exposures



Results


- The UV-C doses required to achieve a 3-log reduction in MRSA and *C. difficile* were 10,000 and 46,000 µJ/cm², respectively
- For both indicators, there was a visible color change from baseline at 10,000 µJ/cm² and a definite final color change by 46,000 µJ/cm² (Figure 1&2)
- Organic load had only a modest impact on UV-C efficacy
- The indicators required only a few seconds to place and were easy to read (Figure 2)

Conclusions

- UV-C doses of 10,000 µJ/cm² and 46,000 µJ/cm² were required to achieve 3 log reductions of MRSA and *C. difficile* spores, respectively.
- The colorimetric indicators provide an easy means to monitor UV-C dosing.
- Additional studies are needed to evaluate use of the indicators in patient rooms including in shaded areas


Acknowledgement

- We thank The Clorox Company and Intellego Technologies for providing devices for testing
- Providing companies did not have any role in planning or design of the study and no funding was received



Results

- The UV-C doses required to achieve a 3-log reduction in MRSA and C. difficile were 10,000 and 46,000 $\mu\text{J}/\text{cm}^2$, respectively
- For both indicators, there was a visible color change from baseline at 10,000 $\mu\text{J}/\text{cm}^2$ and a definite final color change by 46,000 $\mu\text{J}/\text{cm}^2$ (Figure 1&2)
- Organic load had only a modest impact on UV-C efficacy
- The indicators required only a few seconds to place and were easy to read (Figure 2)



Conclusion

- UV-C doses of 10,000 $\mu\text{J}/\text{cm}^2$ and 46,000 $\mu\text{J}/\text{cm}^2$ were required to achieve 3 log reductions of MRSA and *C. difficile* spores, respectively.
- The colorimetric indicators provide an easy means to monitor UV-C dosing.
- Additional studies are needed to evaluate use of the indicators in patient rooms including in shaded areas

Additional hospital studies

- Dr Donskey study presentation – IUVA workshop 14th of January
- *”Ultraviolet-C decontamination of hospital room: Amount of UV light needed”*
- Hospitals already using dosimeters today
- More to come...



Ultraviolet-C decontamination of a hospital room: Amount of UV light needed


Marie Lindblad^{a,d,*}, Eva Tano^b, Claes Lindahl^c, Fredrik Huss^{a,d}

^a Burn Centre, Department of Plastic and Maxillofacial Surgery, Uppsala University Hospital, Sweden

^b Department of Medical Sciences, Section of Clinical Bacteriology, Uppsala University, Sweden

^c Intellego Technologies AB, Sweden

^d Department of Surgical Sciences, Plastic Surgery, Uppsala University, Sweden



INTELLEGO

TECHNOLOGIES

TECHNOLOGY SAVING LIFE AND MONEY