


# Nanoplate Digital PCR for Wastewater Surveillance



**DHS/NIST Workshop: Standards to Support an Enduring  
Capability in Wastewater Surveillance for Public Health**

June 2021

Michael Bussmann, Ph.D.


Associate Director, Global Product Management Digital PCR

## Legal disclaimer


QIAGEN products shown here are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

## Sources of variability in wastewater testing




## Sources of variability in wastewater testing



### Variability in detection

- Process complexity
- Staff skill set
- Partition volume and count

# Detection and absolute quantitation of pathogens in wastewater



## Why digital PCR for wastewater testing?

- A standard curve is not required
- Better inter-laboratory comparability
- Higher precision due to absolute quantification
- Higher robustness for viral detection from complex samples



## Why Nanoplate digital PCR?


- Fully integrated workflow with time-to-result ~2 h
- Easy operating and workflow like qPCR
- Flexible throughput and resolution options
- High multiplexing capability (up to 5plex)
- **Less expertise and experience needed to run digital PCR experiments**

# Accuracy without complexity

Get accurate dPCR results with an easy-to-use and fast system

NIST Reference Material 2372a (human DNA quantification standard)

- Probe based detection (genomic targets are: CTNNB1, EGFR, ERBB2)
- NIST Reference Material was used to generate 10 and 100 copies/uL templates



# Scalable and comparable

## Reliable and reproducible results over Different Systems

- Three runs per instrument and plate types
- Input (expected concentration): 500 copies/μl
- Assay: QIAGEN dPCR Demo Assay (FAM)

<b>Instrument type /Nanoplate 26K 24- well</b>	Mean cp/μl	SD	CV%
QIAcuity One	480	16	3.3%
QIAcuity Four	484	15	3.1%
QIAcuity Eight	480	16	3.4%
Mean	481		
SD	3.83		
CV%	0.80%		

<b>Instrument type /Nanoplate 8.5K 96- well</b>	Mean cp/μl	SD	CV%
QIAcuity One	523	23	4.3%
QIAcuity Four	530	27	4.9%
QIAcuity Eight	521	27	5.1%
Mean	523		
SD	8.86		
CV%	1.69%		

# Scalable and comparable

## Reliable and reproducible results over Different Systems

- Three runs per instrument and plate types
- Input (expected concentration): 500 copies/ $\mu$ l
- Assay: QIAGEN dPCR Demo Assay (FAM)



<b>Instrument type – total</b>	<b>Mean cp/<math>\mu</math>l</b>
QIAcuity One	499
QIAcuity Four	507
QIAcuity Eight	501
Mean	502
SD	22.23
CV%	4.43%



# Predictable partition count

## Robust quantification with high number of valid partitions



- Number of valid partitions per well for Nanoplate 26K 24-well and Nanoplate 8.5K 96-well



# Accuracy in all measurements without being a dPCR expert

## The VPF (Volume Precision Factor)

- Nanoplates and the VPF enable control of the partition volume





- Ensure precision of the quantitation
- Nanoplate batch & reaction/well specific calibration of the partition volume
- Automatically applied

# Example of results: SARS CoV-2 quantification

## Sample: 10-fold serial dilution of positive control ATCC

- Concentration range between  $10^4$  to 1 copy/ $\mu$ l
- Test all samples in triplicates, 4  $\mu$ l template for both
- Panel 1\_3-plex SARS CoV2 N1/N2/E (HEX/Texas Red/FAM)
- “Low/high 10”  $\rightarrow$  10 copies of SARS-CoV-2 target in the presence of  $10^5$  copies each of the fecal indicator targets MS2, HF183 and crAssphage

## N1 (HEX)



## Reference:


Taniuchi, Mami; Associate Professor, Medicine: Infectious Diseases and International Health;  
University of Virginia (Field tester of the QIAcuity One-Step Viral RT-PCR Kit, cus)

# Example of results: Norovirus G1 quantification

## Experimental approach

- Samples – water from swimming pools
- Sample prep – QIAamp Viral RNA Mini kit

Consistent fluorescence intensity in reference channel is a good indicator of uniform partitioning across a nanoplate



Reaction Mix	Target	Sample/NTC/Control	Concentration copies/ $\mu$ L	CI (95%)	Partitions valid	positive	negative	Threshold
B1	Noro g1	23 E2 Janv	0.349	79%	25451	7	25444	60.18
E1	Noro g1	26 E1 Juin	1.9	32.6%	25467	37	25430	110.92
F1	Noro g1	27 E2 Juin	14.0	12%	24651	265	24386	68.85

A background image showing a close-up of a hand reaching out towards the camera, with other hands blurred in the background, suggesting a crowd or a group of people.

Thank you for your attention.

Trademarks: QIAGEN®, Sample to Insight®, QIAcuity®, QIAseq® (QIAGEN Group); Atto™ (ATTO-TEC GmbH); EvaGreen® (BIOTIUM, INC); FAM™, HEX™ (Life Technologies Corporation); Horizon™ (Horizon Discovery Limited); ROX™ (Thermo Fisher Scientific or its subsidiaries); SYBR® (Life Technologies Corporation). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. PROM-18244-001 © 2021 QIAGEN, all rights reserved.