


# Polarization Analyzed SANS

Flip neutron spin at will with respect to the polarization axis using electromagnetic coil of RF flippers.

$^3\text{He}$  cell can spin analyze a divergently scattering beam.



$H$  defines polarization axis.  
Select out one spin orientation.

Structural scattering does not flip neutron spins; magnetic scattering can (provides means to separate the two).

- I. Spin leakage correction (Reduction, J. Appl. Cryst. 45, 546-553, 2012)
- II. Magnetic interpretation of cross-sections (Analysis, J. Appl. Cryst. 45, 554-565, 2012).

## Standard Instrumental Data Corrections

<input checked="" type="checkbox"/> Sensitivity	set DIV file	file: ask
<input checked="" type="checkbox"/> Mask	set MASK file	file: ask
<input checked="" type="checkbox"/> Absolute Scale	set ABS params	parameters: ask

- **2D Detector Sensitivity:** Corrects for bad pixels with DIV file
- **Absolute Scaling:** Normalize to Incident Beam Intensity
- **Mask** scattering from detector edges and “bad” spots

ALSO

### Math for Sample-Specific Corrections

$$\text{COR} = (\text{SAM} - \text{BGD}) - [\text{Tsam}/\text{Temp}](\text{EMP} - \text{BGD})$$

Blocked Beam (BGD) measure  
of electronic and ambient bgr


Sample Scattering files (SAM)

Empty Cell (EMP) measure of bgr  
from sample holder and main beam


Sample and Empty Cell transmissions  
relative to intensity of incident beam

# Polarization Correction

Before correction structural  
Bragg peak dominates





After correction four-fold magnetic  
symmetry becomes apparent



Non spin-flip contains structural scattering and usually is much more intense than magnetic-only spin-flip. Thus, spin-flip is more strongly affected by leakage.

## Getting Started with IGOR for NG7 PASANS Reduction

Igor software and macros found at [http://www.ncnr.nist.gov/programs/sans/data/red\\_anal.html](http://www.ncnr.nist.gov/programs/sans/data/red_anal.html)



### First Steps:

- Pick Path for *GlassyCoFe2O410KHFFData NG7* files
- List files
- Plot 2D data


## Checklist of files needed:

- Open Beam transmission (SM IN, 3HE OUT) for absolute scaling at every detector distance [often obtained during alignment]
- Sample transmission (SM IN, 3HE OUT) needed for absorption scaling
- Blocked beam transmission for each detector distance
- Blocked beam scattering for each detector distance
  
- A Trans. SM OUT, 3He OUT
- B Trans. SM OUT, 3He IN
- } Measures 3He decay (repeat every few hours)
  
- C Trans. (↑↑) SM IN, FLIP OFF, 3HE OFF
- D Trans. (↓↑) SM IN, FLIP ON, 3HE OFF
- E Trans. (↓↓) SM IN, FLIP ON, 3HE ON
- F Trans. (↑↓) SM IN, FLIP OFF, 3HE ON
- } Measures supermirror efficiency, flipper efficiency, and any sample depolarization
  
- G Scatt. (↑↑) SM IN, FLIP OFF, 3HE OFF
- H Scatt. (↓↑) SM IN, FLIP ON, 3HE OFF
- J Scatt. (↓↓) SM IN, FLIP ON, 3HE ON
- K Scatt. (↑↓) SM IN, FLIP OFF, 3HE ON
- } Measures desired four scattering cross-sections

**Measurement conditions:** 1) Sample at 10K and 1.2T with 3.5m and 14m detector settings  
2) Empty at RmT and 0.07T with 3.5m and 14m detector settings



## 1) List files to make two tables



## 2) Associate Sample Trans files with corresponding Open Beam file

TransmissionFiles


T_EMP_Filenames	T_Filenames	T_Labels
JUN16004.SA2_IC	JUN16203.SA2_IC	CD2_SS 0.8 G 300 K Transmission
JUN16004.SA2_IC	JUN16211.SA2_IC	CD2_SS 500 G 300 K Transmission
JUN16021.SA2_IC	JUN16215.SA2_IC	CD2_SS 500 G 300 K Transmission
JUN16022.SA2_IC	JUN16217.SA2_IC	CD2_SS 500 G 300 K Transmission

## 3) Associate Sample Scatt files with corresponding Sample Trans file


ScatteringFiles

S_TRANS_Filename	S_Filenames	S_Labels
JUN16203.SA2_IC	JUN16210.SA2_IC	CD2_SS 0.8 G 300 K Scattering
JUN16211.SA2_IC	JUN16212.SA2_IC	CD2_SS 500 G 300 K Scattering
JUN16211.SA2_IC	JUN16213.SA2_IC	CD2_SS 500 G 300 K Scattering
JUN16211.SA2_IC	JUN16214.SA2_IC	CD2_SS 500 G 300 K Scattering

## Calculate Sample and Empty Transmissions




## 4) Set region of Open Beam to use for normalization. Calculate Transmission.




## Reducing Polarized Data in Four Easy Steps

- 1) Define  $^3\text{He}$  Cells
- 2) Measure  $^3\text{He}$  Time-dependence
- 3) Measure supermirror and flipper efficiencies
- 4) Polarization correct scattering files on absolute scale



## Step 1: Define Your $^3\text{He}$ Cells

- $T_E$  (transmission of glass) is known *a priori*
- $\mu$  (effective cell path length) is known *a priori*
- Nominal values will be available within the IGOR framework.




CellName	lambda	Te	err_Te	mu	err_mu
Maverick	5	0.87	0.01	3.184	0.2
Burgundy	5	0.86	0.01	3.138	0.15
Olaf	7.5	0.87	0.005	3.2	0.005
Hollywood	7.5	0.87	0.005	3.023	0.005
Olaf2	7.5	0.87	0.005	3.2	0.005
OlafB	7.5	0.87	0.005	3.2	0.005
HollywoodB	7.5	0.87	0.005	3.023	0.005
OlafC	7.5	0.87	0.005	3.2	0.005
HollywoodD	7.5	0.87	0.005	3.023	0.005
OlafE	7.5	0.87	0.005	3.2	0.005

**What you need to know:** Every  $^3\text{He}$  analyzer cell has its own decay curve.



## Step 2: $^3\text{He}$ Decay Curve

- $T_{\text{Major}}$  is the transmission of majority spin state
- $\Gamma$  is decay constant in hours
- $P_o$  is initial atomic polarization;  
 $P_{\text{cell}}$  is effective polarization
- May override values, if desired



Use Files A and B: Unpolarized transmissions define time-dependent decay curve.

## Step 3: Polarization Efficiencies

- $P_{\text{Cell}}$  is automatically calculated for times at which polarized transmissions are taken.
- Sample depolarization of the neutron beam is *condition dependent*.

User entered

Field Condition: Cond\_HighField Add Condition ?

R0		Burgundy1									
CondCell_HighF	Cor	Cond_HighField	Cond_HighField	Cond_HighField	Cond_HighField	Cond_HighField	Cond_HighField	Cond_HighField	Cond_HighField	Cond_HighField	Cond_HighField
	UU_Trans	DU_Trans	DD_Trans	UD_Trans	Blocked	PoL_SM_FL	PoL_SM	Include			
Burgundy1	61	62	59	58	63	0.867447	0.885392	1			
Burgundy1	24	25	27	28	63	0.86771	0.889386	1			
Burgundy1	24	25	59	58	63	0.847959	0.900663	0			
Burgundy1	61	62	27	28	63	0.885485	0.872314	0			
Burgundy1	36	37	35	34	63	0.868748	0.88778	1			
Burgundy1	49	50	52	51	63	0.862444	0.88403	1			

Do Average Show Calc Clear Table Clear Row

**AVERAGED RESULTS**

Sam\_depol\*Psm\*Pf 0.866587 +/- 0.018871 Manual Entry

Sam\_depol\*Psm 0.886647 +/- 0.0193037 Snapshot

Use Files C, D, E and F: Measure  $\uparrow\uparrow$ ,  $\downarrow\uparrow$ ,  $\downarrow\downarrow$ , and  $\uparrow\downarrow$  transmissions to obtain polarization of super mirror, polarization of flipper, and sample depolarization.

## Step 4: Polarization Correct Sample and Empty Data

- Associate  $^3\text{He}$  cell used with file
- Can add up to 10 files together
- **Data and empty** scattering files are pol corrected **separately**
- **Fully reduce each detector setting separately**

SAM EMP BGD Condition Cond\_HighField

UU or ++

Run #	Cell
41	Olaf
67	Maverick

DU or - +

Run #	Cell
40	Olaf
68	Maverick
69	Maverick

DD or --

Run #	Cell
30	Olaf
45	Burgundy
90	Maverick

UD or +- -

Run #	Cell
31	Olaf
44	Burgundy
50	Burgundy
53	Burgundy

Load ... Show Coef Matrix Clear Entries

**Pol Correct Data** Change Displayed Data Display 4 X S

Use EMP?  Use BGD?

Sensitivity set DIV file file: Plex\_08Nov2010\_NG3.div,

Mask set MASK file file: DEFAULT.MASK,

Absolute Scale set ABS params parameters: none

Average and Save set AVERAGE params parameters: AVTYPE=Circular;PHI=0;DPHI=0;WID

Reduce Data Save Protocol Recall Protocol

## Followed by Standard Instrumental Data Corrections

**Average and Save:**  
Choose 2D\_Ascii for  
format of corrected data

- Subtraction of Empty Cell
- Detector Sensitivity
- Detector Mask
- Absolute Intensity Scaling

PolCor\_Panel

SAM EMP BGD Condition Cond\_HighField ?

**UU or + +**

Run #	Cell
41	Olaf
67	Maverick

**DU or - -**

Run #	Cell
40	Olaf
68	Maverick
69	Maverick

**DD or - -**

Run #	Cell
30	Olaf
45	Burgundy
90	Maverick

**UD or + +**

Run #	Cell
31	Olaf
44	Burgundy
50	Burgundy
53	Burgundy

Load ...
Show Coef Matrix
Clear Entries

Pol Correct Data
Change Displayed Data
Display 4XS

---

Use EMP?     Use BGD?

Sensitivity        file: Plex\_08Nov2010\_NG3.div,

Mask        file: DEFAULT.MASK,



Absolute Scale        parameters: none

Average and Save        parameters: AVTYPE=Circular,PHI=0,DPHI=0,WID=

Reduce Data
Save Protocol
Recall Protocol

## Final Steps to Generate 1D Data Slices

Generate 1-D Sector Cuts  
(Intensity vs Q)



Splice together data from  
different detector settings

Plot and then **fit with Sasview!**

