



celsis®

RAPID MICROBIAL DETECTION

ABOUT THE SPEAKER



Jonathan Kallay
Senior Director of
Technology and Market
Development

Jonathan (Jon) Kallay is a Senior Technical & Market Development Manager working remotely for the Microbial Solutions product lines. He is a subject matter expert on microbiological investigations for manufacturing facilities that make regulated products. Jon provides practical laboratory experience to help clients identify the optimal path forward for their labs.

Jon received his Bachelor's degree in biochemistry from Denison University before earning a post-graduate diploma in pharmaceutical microbiology from the University of Manchester.



celsis®

RAPID MICROBIAL DETECTION

Sterility Test Flow with Celsis

Sterility Testing



Incubation




Reading Results



Same day Transfer

14 Day Incubation






CELSIS AMPISCREEN® AMPLIFIED ATP BIOLUMINESCENCE

Celsis® ATP-Bioluminescence Reagents

- Two-phase, proprietary enzyme reaction
- All living organisms also contain the enzyme adenylate kinase (AK) as part of their biochemical processes
- Microbial enzymes convert ADP into ATP
- Amplification of ATP levels beyond naturally occurring level.
- Enzymes are not depleted by reaction
- Ability to generate almost unlimited amounts of ATP



the Celsis Accel®
Product Code: 7460288



Complete RMM Solutions.

What does yours look like?

the Celsis Advance II™
Product Code: 7456004



Celsis AMPIScreen®
Pharma Reagent Kit

400 Assay Reagent Kit
Product Code: RST400



Sartorius Sterisart® NF
for Celsis®

Celsis® Qualified Sterisart® NF
Septum with 4 cm dual-needle
Product Code: 16466CR-GSD

Celsis® Qualified Sterisart® NF
Septum with 5.2 cm needle
Product Code: 16467CR-GSD



Hardy Diagnostics Media
for Celsis®

Celsis® Qualified Hardy Diagnostics
Tryptic Soy Broth, USP,
100 mL Septum Top Glass Vial
Product Code: CM1010

Celsis® Qualified Hardy Diagnostics
Fluid Thioglycollate Medium, USP,
100 mL Septum Top Glass Vial
Product Code: CM1015

Charles River Implementation and Validation Support

Celsis® Complete Validation & Reports

Method Validation documentation for
Equivalency, Specificity, Robustness,
Ruggedness, and Limit of Detection.

cGMP Validation Testing in Presence of
Product for Method Suitability, Equivalency,
Limit of Detection, Specificity.

Product Code: VAL6000MF

Celsis® Advantage Reports & Protocols

Method Validation documentation for
Equivalency, Specificity, Robustness,
Ruggedness, and Limit of Detection.


Protocols for Testing in Presence of Product
for Method Suitability, Equivalency, Limit of
Detection, Specificity.

Product Code: VAL6100MF

Installation Support and Ongoing Training

Celsis® 3 Day Initial Installation
and Training
Product Code: TS0003

Celsis® 1 Day Supplemental Training
Product Code: TS0001




Answers the age-old question in
RMMs...

Does it work with CHO?


The Celsis Adapt™ unlocks the ability to use just **one** rapid microbiological method across your entire manufacturing process **from cell culture to final sterility.**

CELSIS ADAPT

Celsis Adapt™ Introduction



Demonstrated detection of broad range of microorganisms, **including fungi.**



Compatible with various cell lines used in bioprocessing of pharmaceuticals and cell & gene therapies.



Aligns with current and new regulatory requirements for the use of alternative rapid microbiological methods.



Sterility Test Flow With Celsis Adapt

Sterility Testing



Incubation Product Cell Removal



Reading Results




Same day Transfer

6 Day Incubation

Same day Testing


CELSIS ADAPT™ PREP PROCESS OVERVIEW

Celsis Adapt™ Methodology & Workflow




Prepare & Incubate

- Prepare sample using compendial method TSB/FTM media used
- Incubate for 3-7 days




Lyse & Concentrate

- Aseptically remove sample aliquot and add to Treatment Solution
- Allow for a minimum dwell time of 15 minutes
- Concentrate using Adapt™



Pipette Concentrate + Celsis LuminASE®

- Combine aliquot of concentrate and reagent to sterile assay cuvettes.




Run on Celsis® instrument

- Load samples into Celsis Advance II™
- Allow 20-minute ATP depletion step
- Obtain rapid qualitative results reported in positive or negative


REDUCTION OF ATP BACKGROUND

Background of 10 mL Cell Sample (CHO, 1×10^6 / mL) + 100 mL TSB/FTM



STEM CELLS METHOD SUITABILITY (PH. EUR. 2.6.27)

1 mL Human Umbilical Stromal Cells (5×10^6 / mL) and Supernatant + 100 mL TSB/FTM, 7 Days Incubation



CELSIS DATA AND VALIDATION SUPPORT

CELL LINE COMPATIBILITY

Celsis Adapt™ Cell

The testing system has undergone rigorous compatibility testing with commonly used cell lines in biopharma processing and cell-based production. This growing list of cell lines include, but are not limited to, the following cell types:

- A9 (Mouse)
- Adipose Stromal Cells (Human)
- ATCC® VR-844™ (Murine Sarcoma Virus)
- BALB 3T3 (Mouse embryo fibroblast)
- CEF (Chicken embryo fibroblast)
- CEL (Chicken embryo liver cells)
- CHO-K1 (Cricetulus griseus, Chinese Hamster ovary)
- CHO-DG44
- CHO-BHK21 (Baby hamster kidney cell)
- HEK-293 (Human embryonic kidney)
- HeLa (Human)
- MDBK (Madin-Darby bovine kidney)
- Mesenchymal stromal cells
- Mrc-5 (Medical Research Council cell strain 5)
- PG4 (S+L- Cat)
- S+L- Mink
- TK6 (Homo sapiens spleen lymphoblast)
- Umbilical Stromal Cells (Human)
- Vero (African green monkey kidney epithelial)
- XC-Rat



UNDERSTANDING THE RELEASE TESTING OPTIONS

Multiple Opportunities to Release Product Faster

3day
RELEASE TEST

Using USP <1071> and Ph. Eur. 2.6.27 approach for ATPMs and short shelf-life products. Utilizes

4day
STERILITY TEST

In accordance with USP <71>, Ph. Eur. 2.6.1, but utilizes a 3rd incubation parameter of TSB at

7day
STERILITY TEST

In accordance with USP <71>, Ph. Eur. 2.6.1, and utilizes compendial incubation parameters.

ALTERNATIVE METHOD VALIDATION REQUIREMENTS

Requirements	Required By	Definition	Charles River Support
POC		Feasibility or principle (i.e., assessing whether the method and accompanying system is suitable for its intended purposes and that it is compatible with the intended product or sample matrix.)	Celsis® Sample Effects and Spiking Studies performed in Applications Lab
Instrument Qualification		Installation Qualification	Celsis Advance II™ Installation Qualification
		Operational Qualification	Celsis Advance II™ Operational Qualification
		Performance Qualification	Celsis Advance II™ Performance Verification Additionally, user must complete internal performance qualification requirements as deemed appropriate by the user
Validation of Alternate Technologies		Specificity	Celsis® Sterility Equivalency (Membrane Filtration) – Report (VAL4000EQ)
		Limit of Detection	
		Equivalence/Comparative Testing/Accuracy	The Celsis AMPiScreen® Rapid Detection Assay Verification of Robustness – Report (TRCELSIS-06.0)
		Robustness	The Celsis AMPiScreen® Rapid Detection Assay Verification of Ruggedness – Report (TRCELSIS-05.0)
		Ruggedness	Dependent upon application; support provided as needed.
		Repeatability	Dependent upon application; support provided as needed.
Method Suitability		Suitability	Method Suitability Test
		Equivalence/Comparative Testing	Equivalence demonstration in presence of product to include relevant environmental isolates

charles river

charles river

charles river

Each customer is responsible for consulting internal quality requirements and Regulatory authority, as appropriate, for validation requirements.

CONTACT US

Jon Kallay

Senior Technical and Market Development Manager
Jonathan.kallay@crl.com

Address:

**251 Ballardvale Street
Wilmington, MA
01887**

Email:

askcharlesriver@crl.com

Website:

www.criver.com

Phone:

877.CRIVER.1