



FORENSICS @ NIST


#NISTForensics

Building the Chemical Foundation for Intelligent Breath Analysis

K. Jeerage, M. Harries, C. Suiter, E. Holland, J. Widegren, T. Lovestead
Fluid Characterization Group
Applied Chemicals & Materials Division, NIST-Boulder

Funding: NIST Special Programs Office

Drug Sampling Options



Blood Sample


1.5 h – 4 h delay



FORENSICS @ NIST

#NISTForensics

Drug Sampling Options



Ethanol

Known correlation between conc. and impairment

Blood-breath ratio \approx 2350 (mean)

Simple elimination profile; not stored in fat tissue

Thermophysical properties – known

Cannabis

Unknown correlation between conc. and impairment


Blood-breath ratio \approx ?

Unknown elimination profile; stored in fat tissue


Thermophysical properties – unknown and difficult to measure



Cannabis is Challenging



ethanol



Δ-9-tetrahydrocannabinol (THC)

metabolites

cannabigerol & cannabinalol


terpenes



FORENSICS @ NIST

#NISTForensics


Cannabis is Challenging



Ethanol Blood Concentration

0 – 3,000,000 ng/mL

SOURCE: Jones et al. *J. Forensic Sci.* **41**, 1996.



THC Blood Concentration

0 – 200+ ng/mL

SOURCE: Huestis *Chem. Biodivers.* **4**, 2007.




FORENSICS @ NIST


#NISTForensics

Properties for Reliable Quantitation

Vapor Pressure



$$P_{\text{SAT}} = \frac{m \cdot R \cdot T}{M \cdot V}$$



Partitioning

FAT ↔ BLOOD ↔ BREATH ↔ SORBENT

$$K_{\text{F/B}} = \frac{[A]_{\text{FAT}}}{[A]_{\text{BLOOD}}}$$


$$K_{\text{B/A}} = \frac{[A]_{\text{BLOOD}}}{[A]_{\text{AIR}}}$$

$$K_{\text{S/A}} = \frac{[A]_{\text{PDMS}}}{[A]_{\text{AIR}}}$$



P_{SAT} – Standard Method

Concatenated Gas Saturation



Time = 1 – 8 weeks

Mononitrotoluene taggants

Widegren et al. *J. Chem. Eng. Data* **55**, 2010.

Low-volatility terpenes

Widegren et al. *Environ. Sci. Technol.* **44**, 2010.

Internal eicosane control

Widegren et al. *Fuel* **90**, 2011

Influence of carrier gas

Widegren et al. *J. Chem. Eng. Data* **60**, 2015.

Anti-oxidants for unstable chemicals

Widegren et al. *J. Chem. Eng. Data* **62**, 2017.





FORENSICS @ NIST

#NISTForensics

P_{SAT} – Rapid Method

Porous Layer Open Tubular
(PLOT) Cryoadsorption



Tetradecane validation

Lovestead et al. *Anal. Chem.* **82**, 2010.

Cannabinoids

Lovestead et al. *Forensic Chem.* **5**, 2017.




FORENSICS @ NIST


#NISTForensics

P_{SAT} – Rapid Method

Porous Layer Open Tubular
(PLOT) Cryoadsorption



Time = 1 – 2 hours



THC: 0.000115 Pa @ 40 °C




FORENSICS @ NIST


#NISTForensics

P_{SAT} – Rapid Method

Porous Layer Open Tubular
(PLOT) Cryoadsorption



Time = 1 – 2 hours



THC: 0.000115 Pa @ 40 °C


Ethanol: 18,000 Pa @ 40 °C
SOURCE: NIST REFPROP Database



FORENSICS @ NIST

#NISTForensics

$K_{\text{Blood}/\text{Air}}$ – PLOT Cryoadsorption Method



Can we capture and concentrate sufficient material
at physiological temperatures?

How does ethanol impact blood/air partitioning?



FORENSICS @ NIST

#NISTForensics

$K_{\text{Fat/Blood}}$ – Non-Invasive NMR Method

Fat Surrogate = Octanol


Blood Surrogate = Water

Eugenol
Signals

Eugenol in D_2O

Octanol + D_2O

Eugenol in
Octanol + D_2O



FORENSICS @ NIST

#NISTForensics


$K_{\text{PDMS/AIR}}$ – Model Feasibility

Is it feasible to predict sorbent-air partitioning with an empirical model?




$$K = \frac{[A]_{\text{PDMS}}}{[A]_{\text{AIR}}}$$

$$n \text{ [moles]} = K * V_s * C_o$$

Partition Coefficient



PDMS Structure



Synthetic Designer Drugs



SPME Fiber



Capillary Microextraction of Volatiles (CMV)




$K_{\text{PDMS}/\text{AIR}}$ – Chemicals for Model

Functional Groups (18)

Non-Ring Groups	Ring Groups	Oxygen Groups
>C<	>C<	-OH (alcohol)
>CH-	>CH-	-OH (phenol)
-CH ₂ -	-CH ₂ -	-O-
-CH ₃		
=C<	=C<	>C=O
=CH-	=CH-	-CHO
=CH ₂		-COO-


Training Set Requirements

- 1) Log K at three temperatures: 60 °C - 180 °C
- 2) Follows van't Hoff Equation with $R^2 > 0.95$



$K_{\text{PDMS/AIR}}$ – Model Evaluation

Compare experimental values (x axis) with predicted values (y axis) for chemicals used to build model (N = 275).



K. Jeerage et al. *20th Symposium on Thermophysical Properties*, June 2018.




FORENSICS @ NIST

#NISTForensics

$K_{\text{PDMS/AIR}}$ – Model Predictions

Predict values for chemicals not used to build model.

Example: C_8 chemicals at 100 °C



ethylidene-cyclohexane

2.112 + 0.002


2.128

0.8 %

Experimental Data

Model Prediction

% Difference



ethenyl-cyclohexane

1.989 + 0.001

1.996

0.4%

K. Jeerage et al. *J. Breath Res.* (in preparation)



FORENSICS @ NIST

#NISTForensics

Summary

Quantitative breath analysis of cannabinoids will be challenging due to their low volatility and their low concentration in blood.

Thermophysical property measurements for large, hydrophobic, low-volatility chemicals push the limits of existing experimental approaches.

Current Partnerships: Prof. Jerry King (University of Arkansas), Prof. Roger Giese (Northeastern University), Prof. April Hill (Metro State University), Prof. Lupita Montoya (University of Colorado)

Outreach: ASTM Committee D37 on Cannabis, Front Range Forensic Chemists Association, Colorado Bureau of Investigation, Cannabis Industry



FORENSICS @ NIST

#NISTForensics

Megan E. Harries Kavita M. Jeerage Tara M. Lovestead

Chris L. Suiter Jason A. Widegren

Boulder Cannabis Research Team

Please contact any of us to discuss working together
or to discuss postdoctoral opportunities.

jeerage@boulder.nist.gov



FORENSICS @ NIST

#NISTForensics