



Nano- and Microtechnology Development for Advanced Scientific Measurements throughout the Solar System

Adrian Southard, Yun Zheng, Stephanie A. Getty

NASA Goddard Space Flight Center

September 13, 2010

Goddard Space Flight Center



NASA Centers (10):

- + Am
- + Dry
- + Gle
- + Go
- + Jet
- + Joh
- + Ker
- + Lar
- + Ma
- + Ste



Nano and Microfabrication Facilities



SWCNT and MWCNT growth



Wet Chemistry Benches

Scanning Electron Microscopy



Electron Beam Lithography



Front/Back-Side Mask Aligners



Low Pressure Chemical Vapor Deposition
September 13, 2010



Physical Vapor Deposition Systems



Reactive Ion Etch Systems

Planetary Targets for *In Situ* Instruments



Sol 20

Sol 24

Mars

Comet Wild 2




Enceladus

Aerogel
Sample
Capture

Mass Spectrometry



- Well established technique for in situ chemical analysis in planetary science missions
- Provides mass-to-charge ratio for every sample constituent – non-specific
- Can be interfaced to complementary analytical techniques for thorough sample characterization, e.g.
 - Pyrolysis (soil)
Thermally evolved volatiles
 - Gas chromatography
(atmosphere, pyrolysis products)
Small, robust molecules
 - Liquid chromatography
(soil extract, liquid)
Complex molecules



Planetary Mass Spec: State of the Art



Gas Chromatograph Quadrupole Mass Spectrometer

- part of Huygens lander
- Cassini mission to Saturn and moons

- Quadrupole mass filter uses AC fields to select transmission at a single mass
- Mass spectrum is acquired by scanning through operating range
- Thermionic filaments are used for electron impact ionization of gas




Planetary Mass Spec: State of the Art




Pyrolysis Gas Chromatograph Quadrupole Mass Spectrometer

- part of Sample Analysis at Mars Instrument Suite
- Mars Science Laboratory rover mission

- QMS similar to Huygens
 - $m \sim 1.3 \text{ kg}$, $P \sim 14.5 \text{ W}$
- Thermionic filaments are used to ionize pyrolysis products or atmospheric gases






In Development: Miniaturized Time-of-Flight Mass Spectrometer



- Time-of-Flight Mass Spectrometer
 - Field emission **Electron Gun** for electron impact ionization
 - Ions are accelerated to a uniform kinetic energy in the **Ion Lens Assembly**
 - Heavier masses travel more slowly than light ones in the **Reflectron** analyzer
 - Arrival of isomass ion packets is registered as a function of time at the **Microchannel Plate Detector**

Ionization for Time-of-Flight Mass Spectrometer





Previous Work: Field Emission in CNTs



- SWCNT, MWCNT, and CNF
 - Various growth techniques


Our approach:
Tower x-section \ll spacing

MWCNTs:
P. G. Collins and A. Zettl (1996)




SWCNTs:
J.-M. Bonard, *et al.* (1998)

CNFs:
E. Minoux *et al.* (2005)



Patterned MWCNTs:
S. Fan *et al.* (1999)




Selected References

- A. G. Rinzler *et al.*, *Science* **269**, 1550 (1995).
- W. A. de Heer, A. Chatelain, D. Ugarte, *Science* **270**, 1179 (1995).
- P. G. Collins and A. Zettl, *Appl. Phys. Lett.* **69**, 1969 (1996).
- Q. H. Wang *et al.*, *Appl. Phys. Lett.* **70**, 3308 (1997).
- J.-M. Bonard, *et al.*, *Appl. Phys. Lett.* **73**, 918 (1998).
- W. Zhu *et al.*, *Appl. Phys. Lett.* **75**, 873 (1999).
- S. Fan *et al.*, *Science* **283**, 512 (1999).
- X. Xu and G. R. Brandes, *Appl. Phys. Lett.* **74**, 2549 (1999).
- N. de Jonge and J.-M. Bonard, *Phil. Trans. R. Soc. London A* **362**, 2239 (2004).
- E. Minoux *et al.*, *Nano Lett.* **5**, 2135 (2005).

Previous Work: Field Emission in CNTs



W. A. de Heer, A. Chatelain, D. Ugarte (1995)




- Fowler-Nordheim Tunneling
 - Field enhancement factor
 - Reported values 400-1200 (for MWCNTs)

$$J = K_1 E^2 \exp\left(-\frac{K_2}{E}\right)$$



$$K_2 = \frac{B\Phi^{3/2}}{\beta}$$

Carbon Nanotube Electron Gun



Carbon Nanotube field emitters


- Low power (100x < thermionic)
- Scalable for high sensitivity, redundancy
- Efficient ionization of the parent species
- Challenges: Lifetime and stability



Recent Work: Assembly and Integration




Cathode-grid integration



CNT E-gun for Mini Mass Spec



- Low operating voltages achieved for efficient gas ionization
- Fowler-Nordheim behavior confirmed
 - Field enhancement factor ~ 900
- Microamps of transmitted current in triode mode
- Persists for several hundred hours in high vacuum




CNT E-gun for Mini Mass Spec



MEMS Integration of CNT e-gun for clean packaging


- Long lifetime
- Reduced current noise
- Reliable fabrication



September 13, 2010

A. Southard/NASA GSFC



TOF-MS Integrated Testing



Lens Optimization




Focal plane coincident with
ionization region



Integrated TOF-MS Performance



- Mass resolution > 270




Future Steps




- Opportunities for TOF performance improvement:
 - **Increased sensitivity**
 - Electron current* –
 - Better transmission through revised geometry
 - Optimized emitter geometry
 - Pulsed e-gun
 - Ion transmission* –
 - Improve field uniformity to increase ion yield
 - **Increased mass resolution**
 - Reduce ionization volume
 - Minimize rise time on pulse electronics
 - **Improved emitter lifetime**

Liquid Chromatography-Mass Spectrometry



Miniaturized liquid chromatograph





Successful lab techniques → planetary surface

- Nanoelectronic charge displacement detector at output
- Chemicals elute at different rates based on interactions with LC column
- Data = retention time spectrum

Chemical Separation by Liquid Chromatography





Laser-induced Fluorescence



Design of Micro-LC column



- Packing of microbeads to form stationary phase
- Longer channels → can better distinguish retention times, but higher pressure requirements
- Designed micro channels with varying lengths, 40 mm-100 mm using wafer scale processing.





ChemFET as charge displacement detector



Advantages:

- Fully electronic
 - Microwatts
- Non-specific
 - ChemFET can detect any organic species...
 - ...but requires careful calibration
- Nanoscale
 - Redundancy can reduce risk
- Integratable with Si-based technologies
 - Can enable new on-chip functionality (measure pH, flow rate, buffer concentration)




Electronic Means of Chem/Bio Detection



Requirements

- Sensitive/Selective
- Low Power and Compact
- Electronically Addressable
- Autonomously Operable
- Robust and Reliable






Parameter	SiNW
<i>Metallic or Semiconducting</i>	Semiconducting
<i>Sensor-DNA probe</i>	Direct surface attachment
<i>Channel orientation</i>	Aligned on array
<i>Growth temperature</i>	450°C

- Silicon Nanowires*
- Nanoscale diameter
 - Single-NW limit attainable
 - Scalable for manufacturability and redundancy

SiNW ChemFET




- SiNW fabrication – wafer-scale manufacturability demonstrated
 - Wafer-scale pillar growth for control of nanowire placement
 - Near 100% NW device yield
 - Excellent device uniformity



Transistor = current modulation with gate voltage

Electrospray ionization




http://www.magnet.fsu.edu/education/tutorials/tools/ionization_esi.html

Addition of a lens to the ESI nozzle




Design features:

- Can be tested with a commercial mass spectrometer and is compatible with microLC and TOF MS.
- Utilizes micro-nozzle to enhance atomization of jet.
- Electrodes facilitate formation of an ion beam
- Achieves high electric field for efficient charging of analyte without using high voltages



New design of ESI nozzle



Future Work



- Explore SiNW surface functionalization
 - Specificity
 - Multiplexing
- Liquid sample extraction
 - Extraction from solid samples
 - Concentration of analyte
- Miniaturized liquid sample handling components e.g. valves, pumps
- Higher aspect ratio nozzles



Highly interdisciplinary work
End-to-end systems development

- But this is just one facet of technology development at NASA...
 - Earth-observing instruments
 - Astronomical instruments
 - Spacecraft structures
 - Propulsion
 - Human spaceflight

Acknowledgements



GSFC Astrochemistry Lab

- Jason Dworkin
- Mildred Martin

GSFC Atmospheric Experiments Lab

- Paul Mahaffy
- Daniel Glavin
- Jennifer Eigenbrode
- William Brinckerhoff

GSFC Planetary Magnetospheres Lab

- Bill Farrell
- Telana Jackson
- Haydee Aguilar

CUA/Astrochemistry Lab

- Gunther Kletetschka
- Tomoko Adachi
- Vilem Mikula

GSFC Materials Engineering Branch

- Todd King
- Dan Stewart
- Dewey Dove
- Bruno Munoz
- Rachael Bis
- Jonathon Brame (Summer Intern)
- Nathan Woods (Summer Intern)

GSFC Detector Systems Branch

- Yun Zheng
- Manuel Balvin
- Carolina Gallardo
- David Franz
- Harold Isenberg
- Carl Kotecki
- Mary Li
- Nick Costen
- Larry Hess
- Patrick Roman

UMCP

- Mel Gomez
- Konrad Aschenbach

NCI ATC

- Javed Khan
- Belhu Metaferia
- Jun Wei
- Young Song

NIH/NCI

- Piotr Grodzinski
- Eric Hale
- Jerry Lee
- Nicholas Panaro

This work was supported by the NCI NanoBioSensor Initiative and the GSFC Internal Research and Development Program



Please contact Stephanie Getty for more information or opportunities to collaborate.

**Stephanie Getty
Atmospheric Experiments Laboratory
Materials Engineering Branch
Stephanie.A.Getty@nasa.gov
301-286-9760**

Thank you for your attention.