

# The use of microfluidic and droplet-based digital PCR platforms for DNA quantitation

Standardisation of Genome Amplification Technologies Meeting  
29 May 2014

**Jo Lynne Harenza, Ph.D.**  
Applied Genetics Group  
Biomolecular Measurement Division

**NIST**

**National Institute of  
Standards and Technology**  
U.S. Department of Commerce

# Background on Digital PCR

*Proc. Natl. Acad. Sci. USA*  
Vol. 96, pp. 9236–9241, August 1999  
Genetics


1st paper on digital PCR

## Digital PCR

BERT VOGELSTEIN\* AND KENNETH W. KINZLER


The Howard Hughes Medical Institute and the Johns Hopkins Oncology Center, Baltimore, MD 21231

- Diluted DNA is partitioned into many volumes, ideally containing 0 or 1 template, thermalcycled, and **accessible** targets are counted




# Background on Digital PCR

- Poisson statistics are used to determine an starting DNA copy number without the use of a calibrant
- Derived units, **counts per  $\mu\text{L}$** , are traceable to the SI unit, count one and the meter ( $1 \text{ L} = 10^{-3} \text{ m}^3$ )



# Background on Digital PCR

- Poisson statistics are used to determine an starting DNA copy number without the use of a calibrant
- Derived units, **counts per  $\mu\text{L}$** , are traceable to the SI unit, count one and the meter ( $1 \text{ L} = 10^{-3} \text{ m}^3$ )



**Linearity between ~10-50 % positive droplets/chambers increases accuracy of measurement by decreasing the standard error**




# Droplet digital PCR (ddPCR)

## QX100/200 by Bio-Rad Laboratories

- 96 well PCR plate; 10,000-20,000 droplets per sample
  - Bio-Rad estimates: {0.89, 0.91, 1.0} nL droplets
  - NMI Australia: 0.868 +/- 2 % nL droplets
- end-point DNA quantitation

### Fill Cartridge & Generate Droplets



### Thermal cycle




### Read Droplets



# ddPCR Data Output

samples

NTCs




# Droplet digital PCR (ddPCR)

## RainDrop™ System by RainDance Technologies

- 8 panel chip
- up to 1 million droplets per sample — advantageous for detecting rare mutations
- real-time data collection and droplet size estimation during droplet generation

**Fill chip with samples**



**Generate Droplets**



**Thermal cycle & Read Droplets**




# Microfluidics: chamber digital PCR (cdPCR)

## Biomark™ HD System by Fluidigm Corporation

- 12x765 (6 nL) or 48x770 (0.85 nL)
  - 9,180 or 36,960 chambers
- Real-time data collection at every PCR cycle

**Fill chip with samples**



**Load samples into chambers**




**Thermal cycle**



# cdPCR Data Output & Analysis


chamber heat maps



sample


NTC

real-time PCR curves



ogives

(frequency distribution plots)




(Jo Lynne Harenza, Dave Duewer)



# cdPCR Data Output & Analysis


chamber heat maps



sample


NTC

real-time PCR curves



ogives

(frequency distribution plots)



(Jo Lynne Harenza, Dave Duewer)


# Interesting Findings

- Artifacts on the Bio-Rad ddPCR system
- Assay Optimization: qPCR v. digital PCR
- Plasmid DNA: digested v. supercoiled plasmid DNA
  - Attempt to relax supercoiled plasmid DNA with DMSO
- Assays with multiple targets are revealed with dPCR
- Master mix contamination




# Artifacts of ddPCR

- **Shifting** of positive and negative droplets
  - Larger than average size droplets generated - caused by a problem with the first generation of cartridges
  - Results in concentration inflation — omit samples from analysis





# Artifacts of ddPCR

- Droplet **shearing**
  - Happens during pipetting of droplets — either into the cartridge or into the PCR plate
  - Can occur when bubbles in the sample mixture rise within the tip and burst/shear droplets
  - affect?




# Temperature gradients are still required to optimize dPCR assays

- For ddPCR, want separation of positive and negative droplets
- For cdPCR, want efficient amplification at a reasonable  $C_T$  and good curve morphologies





(Erica Butts & Margaret Kline)

# Efficient qPCR assays are not always good assays for dPCR






## human *ACTB*

### cdPCR



# Linear v. Supercoiled DNA


- Supercoiled DNA contains many “**late-starters**” and therefore, the target never fully amplifies





# DMSO relaxation of supercoiled pDNA




- DMSO is commonly used to relax plasmid DNA in qPCR [3, 5, 7, 9%]
- DMSO reduced the efficiency of cdPCR reactions
  - Concentrations of 7% and 9% completely inhibited amplification



(Jo Lynne Harenza, Dave Duewer)

# Assays with multiple targets can be detected with digital PCR

## human *RPLP0*



two distinct populations being amplified

two distinct populations being amplified

(Erica Butts & Margaret Kline)



# Many commercial master mixes contain leftover plasmid DNA from recombinant *Taq* production

JOURNAL OF CLINICAL MICROBIOLOGY, Jan. 2005, p. 530–531  
0095-1137/05/\$08.00+0 doi:10.1128/JCM.43.1.530–531.2005  
Copyright © 2005, American Society for Microbiology. All Rights Reserved.

Vol. 43, No. 1

## Presence of $\beta$ -Lactamase Gene TEM-1 DNA Sequence in Commercial *Taq* DNA Polymerase

Biotechnology Letters (2006) 28: 321–325  
DOI 10.1007/s10529-005-5931-3

© Springer 2006

---

### The occurrence of antibiotic resistance genes in *Taq* polymerases and a decontamination method applied to the detection of genetically modified crops


André Perron, Philippe Raymond\* & Robin Simard  
*St-Hyacinthe Laboratory, Canadian Food Inspection Agency, Casavant Blvd West, J2S 8E3 3400,  
St-Hyacinthe, Quebec, Canada*

# Negative Control Contamination


qPCR Taqman

**Assay 1**

cdPCR





NTC  
amplification



**Assay 2**

No  
NTC  
amplification



# Master mixes can be decontaminated

- ArcticZyme PCR Decontamination Kit — DNase specific to dsDNA, followed by heat inactivation with DTT
- Restriction enzyme digestion of dsDNA such that the target is no longer intact (requires *a priori* knowledge of contaminant sequences)<sup>#</sup>
- UV irradiation — sequence and length-dependent<sup>\*</sup>
- $\gamma$ -irradiation induces dsDNA breaks, but high levels destroy *Taq*'s enzymatic activity<sup>^</sup>


<sup>#</sup>Sharma, et al, 1992. A simple method for elimination of unspecific amplifications in polymerase chain reaction. *Nucleic Acids Res.* **20**: 6117–6118.

<sup>\*</sup>Corless, et al, 2000. Contamination and sensitivity issues with a real-time universal 16S rRNA PCR. *J. Clin. Microbiol.* **38**: 1747–1752.

<sup>^</sup>Champlot, et al, 2010. An Efficient Multistrategy DNA Decontamination Procedure of PCR Reagents for Hypersensitive PCR Applications. *PLoS ONE*, **5**(9), e13042.

# Negative control correction is also valid for contaminated master mixes


- A linearity study can be used to validate NTC correction
- Before correction, the regression line is appreciably curved
  - Compared an assay free of NTC contamination with one exhibiting contamination



(Jo Lynne Harenza, Dave Duewer)

# Negative control correction is also valid for contaminated master mixes


- A linearity study can be used to validate NTC correction
- Before correction, the regression line is appreciably curved
  - Compared an assay free of NTC contamination with one exhibiting contamination



(Jo Lynne Harenza, Dave Duewer)


# Negative control correction is also valid for contaminated master mixes

- A linearity study can be used to validate NTC correction
- Before correction, the regression line is appreciably curved
  - Compared an assay free of NTC contamination with one exhibiting contamination



(Jo Lynne Harenza, Dave Duewer)

Use multiple assays and treating results interlaboratory study to determine at consensus value



(Margaret Kline, Ross Haynes, Jo Lynne Harenza, Dave Duewer)



# Summary

- Useful information can be gained from combining cdPCR with ddPCR
  - It is important to examine droplet and chamber traces to uncover artifacts or amplification biases
- Concordance between digital PCR instruments can be achieved following proper assay optimization
- Digital PCR measures only **accessible** targets and therefore, may underestimate the true DNA quantity in certain samples
- Digital PCR quantitation is SI-traceable

# Acknowledgements

## **Applied Genetics Group - Clinical Effort**

Peter Vallone, Ph.D. - Group Leader

Jo Lynne Harenza, Ph.D.

Ross Haynes, M.S.

Margaret Kline, M.S.

## **Chemical Sciences Division - Statistics**

David Duewer, Ph.D.

