



# Yeast Cells as a Candidate Reference Material to Support Training in On-Site Biological Agent Sampling and Detection

Nancy J. Lin, Sandra M. Da Silva, James J. Filliben  
August 14, 2017

*for the NIST-DHS-FDA Standards for Pathogen Detection Workshop*

# Acknowledgement and Disclaimers

The Department of Homeland Security (DHS) Science and Technology Directorate funded the production of this material under Interagency Agreement FTST-17-00017 with the National Institute of Standards and Technology (NIST).

## *Disclaimers:*


- *All opinions expressed in this presentation are the author's and do not necessarily reflect the policies and views of DHS, NIST or affiliated venues.*
- *Certain commercial equipment, instruments, or materials are identified in this paper in order to specify the experimental procedure adequately. Such identification is not intended to imply recommendations or endorsement by NIST, nor is it intended to imply that the materials or equipment identified are necessarily the best available for the purpose.*

# “The biological threat is real and growing...”



A National Blueprint for Biodefense, Oct. 2015  
Bipartisan Report of the Blue Ribbon Study Panel on Biodefense

## *Malicious intent*



## *Benign human activity or simple change of nature*



# Standards to Support Field Biological Agent Detection

*DHS S&T – NIST Interagency Agreement*


**Goal:** Develop standards and methods to support field biothreat detection and biosurveillance

- **Microbial reference materials** for training
- Methods, metrics and standards to characterize biological test materials
- Documentary standards to support field response mission capability

**Impact:** Increasing confidence in field results and improving National ability to detect and respond to suspected biological incidents



*Components of a Biothreat Field Response Capability*



# Standards to Support Field Biological Agent Detection

DHS S&T – NIST Interagency Agreement

**Goal:** Develop standards and methods to support field biothreat detection and biosurveillance

- Microbial reference materials for training
- Methods, metrics and standards to characterize biological test materials
- Documentary standards to support field response mission capability

**Impact:** Increasing confidence in field results and improving National ability to detect and respond to suspected biological incidents



Poster by Nate Olson


“Challenging a bioinformatic tool’s ability to detect microbial contaminants using *in silico* whole genome sequencing data.”

# Most Training Uses Biothreat Agents or Near Neighbor Organisms

- Health and safety risks
- Need for specialized facilities
- Limited material availability
- False positives during real events from contamination
- False positives during training from the environment


# Surrogates Needed: Non-threat, biological materials

Evaluate, challenge, and establish confidence in biological assessment in the field



# Format of a Surrogate Material

## Example workflow



## Desirable Properties


- Well-characterized, quantified
- High concentration ( $>10^7$  cells)
- Long (multi-year) shelf-life
- Stable at 25 °C or 4 °C
- Amenable to powder formation, aerosolization, etc.
- Low cost to end users
- Inactivated?



# Objective



Develop, characterize, and demonstrate modified yeast cells as a surrogate for biothreat agents

- Stable, versatile whole cell material
- Quantification of total cells for DNA-based detection
- Relevant protocols that enable users to expect positive detection



# Develop: Modified *Saccharomyces cerevisiae* NE095

- Designed to challenge nucleic-acid based detection technologies
- Procured lyophilized yeast pellets from Microbiologics, Inc.
  - Verified no PCR inhibition due to lyophilization matrix



\*External RNA Controls Consortium (ERCC)

- *External RNA Controls Consortium (ERCC) DNA sequences are part of NIST Standard Reference Material (SRM) 2374: DNA Sequence Library for External RNA Controls.*
- *ERCC-00095 corresponds to the latter three (of eight total) open reading frames in the phosphate specific transport complex component of *M. jannaschii*.*

# Lyophilized *S. cerevisiae* NE095

Minimized real and perceived risk




Readily available, can use almost anywhere



No false positives in real events from equipment contamination




Indicator of a broken process



Low DNA extraction efficiency to challenge the process




DNA target eliminates environmental false positives



Lyophilized yeast can be crushed into a powder



Quantitative material can track performance




# Initial Studies Demonstrated Successful Detection by Potential End Users


- Interlaboratory study with potential end users detected yeast cells in suspension\*
- 4<sup>th</sup> CST confirmed DNA extraction from the yeast via multiple methods

## \*Interlaboratory study participants

- Florida Dept. of Health
- Michigan Dept. of Community Health
- Minnesota Dept. of Health
- New York Department of Health
- Washington Dept. of Health
- 4<sup>th</sup> Civil Support Team (CST), Georgia Army National Guard



Da Silva et al, *Biomol Detect Quant*, 2016



# Characterize: Quantify Yeast Cells

Quantity (total cells)

**Coulter counter**

Hemocytometer



**Reference value**

Viability (live cells)

**Plate counting**

Live/dead staining

DNA insert stability

qPCR



WGS

- Homogeneity
- Stability
- Fitness for Purpose

# Total Cells per Vial


## Coulter Technology

## Imaging-Based



No significant difference in the mean values  
from the two methods.

# Significantly Higher Standard Deviation for Imaging-Based Method




# Characterization of Prototype Batch

## Homogeneity

Measurement	Cells per vial x 10 <sup>7</sup>	Vials	Reps
Total cells (Coulter)	3.81 ± 0.51 (13.3 %)	28	2
Viable cells (Plating)	0.095 ± 0.018 (18.9 %)	14	1





Viability =  
2.50 ± 0.58 %

## Real-time stability (>2 years)



Total and viable  
cell number  
unchanged

## Accelerated stability (~4 months)



Total cell number  
unchanged

At 50 °C, viability  
drops to ≈ 0



# Demonstrate: Operation Vigilant Sample IV – July 2015

- Sample collection and biological detection exercise conducted in real-time
- Designed to help define a national exercise template for NGB CST Commanders
- Obtained EPA approval to use the yeast
- Led by CPT Bryon Marsh, 4<sup>th</sup> CST
  - 4<sup>th</sup> and 48<sup>th</sup> NGB CSTs
  - FL, GA Dept. of Health (CDC LRN)
  - BioWatch (DHS OHA)
  - EPA
  - Local responders



# Yeast Incorporation into the Exercise


Day 1

On-site yeast prep by CSTs



Day 2

Team assembly and JIT training




Plume model



Day 3

Sample Placement



**Yeast successfully detected.**

Sample collection




Decon and chain of custody



On-site DNA extraction and qPCR



On-site results



# Ongoing Activities

## Yeast as a surrogate powder

- LOD study with PHLs - data being analyzed
- Interlab study on sampling and detection of yeast powder with first responders and LRN (upcoming)



## Yeast as a NIST RM

- Reference value based on total (not viable) cells
- EPA Microbial Commercial Activity Notice required



## ASTM standards

- WK42642 guidance on surrogate materials
- Draft standard for specific applications



# Summary

- Demonstrated that yeast material can be used in place of biothreat agents for training and workflow assessments in the field
- Validated a method to quantify total yeast cells (Coulter)
- Identified a yeast format that is stable and versatile (lyophilized)
- Developed a robust protocol for field training using yeast dried onto a surface
- Paving the way for a first-of-its-kind NIST RM
- The yeast material is a critical part of the developing Quality Assurance infrastructure to support reliable, consistent results from the First Responder Community

## Broader Applicability

- Need for whole-cell based reference materials
- Lessons learned from the yeast apply to other microorganisms, with some caveats
- Next steps:
  - Validate methods to quantify total bacterial cells
  - Develop reference materials based on whole bacterial cells?
  - Extend the mixed microbial gDNA candidate RM for biosurveillance applications




Sandra Da Silva



Jim Filliben

Supported by  
Interagency Agreement  
FTST-17-00017 between  
DHS S&T and NIST



Jayne Morrow



Nate Olson



Lindsay Harris

Steve Lund (NIST)

Zvi Kelman (NIST)

Capt. Bryon Marsh (4<sup>th</sup> CST, ret.)

SFC David Rygmyr (4<sup>th</sup> CST)

Jeremy Clancy (Howard Co. Fire Dept.)

Association of Public Health Laboratories (APHL)

Chris Mangal (APHL)

CDC Laboratory Response Network (LRN)

Rich Ozanich (PNNL)

Yonas Nebiyeloul-Kifle (DHS S&T)