



# Small angle scattering study of NIST SRM1549a: **Bovine Whole Milk Powder**

**Lisa Hong (U. Maryland Baltimore County)**

Summer Undergraduate Research Fellowship  
at the NIST Center for Neutron Research




Image from: <http://www.borrisoleigh.ie/2009-11-18/>

## To contribute towards:

- a better assignment of the casein micelle scattering profile

→ a better understanding of structure and applications ←

Typical Small Angle Scattering of Milk



Adapted from Ingham et. al (2016). *Soft Matter*. 12: 6937-6953

# Milk has been part of the human diet for ~8000 years. What's new?

**We don't fully understand it yet**


- None of the current models explains all experimental results for the particles that stabilize milk : **CASEIN MICELLES**
- NIST SRM 1549a was the chosen sample for this project:
  - certified mass fraction values
  - contributions from this project's findings






# What are casein micelles?

Sub-micron particle with various sizes depending on many factors



## Calcium Phosphate Nanoclusters (CPN)



# Function/Applications of Casein Micelles


## Biom mineralization

- Milk is supersaturated in calcium phosphate. Caseins stabilize it.

## Drug & Nutrient Carriers

- Biocompatible, biodegradable

## Inhibit amyloid fibrils



## Dairy quality & performance

- Micelle size influences processing and texture of dairy products

# Research Strategy

**Stabilize**

**Simulated Milk Ultra Filtrate**



- Buffer mimics milk
- Added natural stabilizers
- Temperature control


**Simplify**

Skim the milk


Fractionate the Casein Micelles

Use SANS w/  
contrast variation:

● CPN



CM












# Simulation Fuzzy Sphere R=150 nm (CM)

# Why neutrons?

# Simulation CPN Nanoclusters R= 4 nm



Monodisperse


Polydisperse

Monodisperse


Polydisperse




# Why neutrons?



$\text{H}_2\text{O}$



70%  $\text{D}_2\text{O}$




98%  $\text{D}_2\text{O}$



# Small Angle Neutron Scattering (SANS)

SANS is the least intrusive, most informative technique

- Covers broad range of angles
- Does not require cooling, additives, etc. (less artifacts from sample preparation)



UV-Visible Signal easily saturated

Used instead:


- UV-vis of diluted fractions
- Dynamic Light Scattering

# Fractionation helps – Optimization of SEC runs


Faster flow, 2ml fractions, more fractions pooled together. Run at 4°C

Less column saturation, 0.5 ml fractions, less volume pooled together. Run at 10°C

## 1<sup>st</sup> COLUMN





## 2<sup>nd</sup> COLUMN



SANS profiles in  $D_2O$  buffers have better signal-to-noise but also match out the CPN contributions. The change in size distributions is however still clear.

# Fractionation helps – Contrast Variation






# Conclusions

- SEC can reduce the natural polydispersity of milk for improved resolution
- Mid-q feature of the spectrum is matched out in D<sub>2</sub>O buffers (CPN matched out).
- USANS-SANS provides information on the different components but sample monodispersity can only be kept for 1-2 days.

**\*VSANS please!\***





# Acknowledgments

- S. Teixeira (SURF Mentor, U. Delaware/NIST)
- T. Hoopes (IBBR), T. Cleveland, & M. Shremshock (NIST)
- L. Wood (MML Material Measurements Laboratory)
- B. Brooks (RET program)
- NCNR SURF team: J. Dura & J. Borchers






Image from: <https://www.snorgtees.com/milk-i-am-your-father>



## DISCLAIMER

Certain commercial equipment, instruments, or materials were identified to foster understanding. Such identification does not imply recommendation or endorsement by the National Institute of Standards and Technology, nor does it imply that the materials or equipment identified are necessarily the best available for the purpose.









Image credit: Susana Teixeira (NIST)






Peter Hristov et. al (2016). Measurement of Casein Micelle Size in Raw Dairy Cattle Milk by Dynamic Light Scattering, Milk Proteins - From Structure to Biological Properties and Health Aspects

Holt et. al (2012). Darwinian transformation of a 'scarcely nutritious fluid' into milk. *Journal of Evolutionary Biology*









Peter Hristov, et. al. (2016). Measurement of Casein Micelle Size in Raw Dairy Cattle Milk by Dynamic Light Scattering, Milk Proteins - From Structure to Biological Properties and Health Aspects

## •Size Exclusion Chromatography: Separation by Size

- Akta Purifier System
- HiPrep 26/60 Sephacryl S-500 Column (Matrix: cross-linked copolymer of allyl dextran and N,N'-methylene bisacrylamide)

