

## Specifications

	Evolution 201 UV-Visible Spectrophotometer	Evolution 220 UV-Visible Spectrophotometer
Optical Design	Double-beam with sample and reference cuvette positions; Czerny-Turner Monochromator	Double-beam with sample and reference cuvette positions; Application Focused Beam Geometry; Czerny-Turner Monochromator
Spectral Bandwidth(s)	1.0 nm	Variable: 1 nm; 2 nm; AFBG Microcell optimized; AFBG Fiber optic optimized; AFBG Materials optimized
Light Source	Xenon flash lamp, 3-year warranty (5 years typical lifetime)	
Detector	Dual Silicon Photodiodes	
Scan Ordinate Modes	Absorbance, % Transmittance, % Reflectance, Kubelka-Munk, log (1/R), log (Abs), Abs*Factor, Intensity	
Resolution	>1.6 (peak-to-valley ratio; toluene in hexane)	
Wavelength		
Range	190–1100 nm	
Accuracy	±0.8 nm (full range 190 to 1100 nm) ±0.5 nm (546.11 nm mercury line)	
Repeatability	≤0.1 nm (546.11 nm mercury line, SD of 10 measurements)	
Scanning Speed	<1 to 6000 nm/min; variable	
Data Intervals	10, 5, 2, 1.0, 0.5, 0.2, 0.1 nm	
Photometric		
Range	>3.5 A	
Display Range	-0.3 to 4.0 A	
Accuracy – Instrument	0.5 A: ±0.004 A 1A: ±0.006 A 2A: ±0.010 A Measured at 440 nm using neutral density filters traceable to NIST/NPL	
Accuracy – Sealed Solutions (EP/BP/TGA)	±0.010 A (60 mg/L K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> )	
Noise	0A: ≤0.00015 A 1A: ≤0.00050 A 2A: ≤0.00080 A 260 nm, 1.0 nm SBW, RMS	
Drift (Stability)	<0.0005 A/hr 500 nm, 1.0 nm SBW, 1 hour warm-up	
Stray Light	KCl, 198 nm: ≤1% T NaI, 220 nm: ≤0.05% T NaNO <sub>2</sub> , 340 nm: <0.05% T	
Baseline Flatness	±0.0010 A 200–800 nm, 1.0 nm SBW, smoothing	
Keypad	Sealed Membrane	
Local Control Option		
Display	Touchscreen LCD panel; 800 × 480; 17.8 cm (7 in) diagonal	
Operating System	Microsoft Windows XP embedded	
Dimensions	62.2 cm L × 48.6 cm W × 27.9 cm H (24" L × 19" W × 11" H)	
Weight	14.4 kg (32 lb)	
Electrical Supply	100–240 V, 50–60 Hz, selected automatically 150 W maximum	

[www.thermoscientific.com/uv-vis](http://www.thermoscientific.com/uv-vis)

©2010, 2011 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. Windows is a registered trademark of Microsoft Corporation. Spectralon is a registered trademark of Labsphere, Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Electron Scientific Instruments  
LLC, Madison, WI USA is ISO Certified.

**Africa-Other** +27 11 570 1840  
**Australia** +61 3 9757 4300  
**Austria** +43 1 333 50 34 0  
**Belgium** +32 53 73 42 41  
**Canada** +1 800 530 8447  
**China** +86 10 8419 3588  
**Denmark** +45 70 23 62 60

**Europe-Other** +43 1 333 50 34 0  
**Finland/Norway/Sweden**  
+46 8 556 468 00  
**France** +33 1 60 92 48 00  
**Germany** +49 6103 408 1014  
**India** +91 22 6742 9434  
**Italy** +39 02 950 591

**Japan** +81 45 453 9100  
**Latin America** +1 561 688 8700  
**Middle East** +43 1 333 50 34 0  
**Netherlands** +31 76 579 55 55  
**New Zealand** +64 9 980 6700  
**Russia/CIS** +43 1 333 50 34 0  
**South Africa** +27 11 570 1840

**Spain** +34 914 845 965  
**Switzerland** +41 61 716 77 00  
**UK** +44 1442 233555  
**USA** +1 800 532 4752

BR51944\_E 12/11M

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific