


Visualizing Insulin Structure Under Shear Stress

Andre Chen

Mentor: Grethe Jensen

NIST Center for Neutron Research



Insulin is a hormone that **regulates glucose intake** by cells.

Patients with diabetes are **insulin-deficient**.


Diabetes is treated by injection of **engineered insulin analogues**.

INSULIN THERAPEUTICS



OBJECTIVE

Long-acting insulin drug for once-daily dosing


Hexamers



Rod-like Structures



Next Step: Does Shear Change Insulin Structure?



Subcutaneous Injection


Shear rates on order of 10^5 - 10^6 s^{-1}

Blood Vessels


Shear stresses as high as 9.8 N/m^2

Possible Outcomes



Rod Lengthening/Breaking



Alignment with Applied Shear



Interactions



Project Overview

Objective: Determine if shear affects insulin self-assembly



Insulin samples prepared:

Insulin concentrations (mM): **0.6** | 1.2 | 1.8

NaCl concentrations (mM): **30** | **150** | 450

Formulation

Physiological

Screens
electrostatic
repulsions




Samples sheared at different shear rates to probe
their **viscosity**



Small-angle neutron scattering applied during
shear for structural data

Viewing Nanostructure With Small Angle Neutron Scattering




Why SANS?


Probes appropriate length scale to view:

- Structure
- Orientation of structures



5 Hexamers per rod




10 Hexamers per rod



20 Hexamers per rod



Isotropic orientation



Alignment in flow






Image from <http://cns.che.udel.edu/>

Rheology Measurements



Instrument: Anton Paar
MCR 502 Rheometer


Shear Rates: 0-3500 s^{-1}



**Titanium
Couette Cell**

**Sample
Solution**

Front View




**Tangential
Velocity
Profile**

Top View

Images:
Rheometer:
<http://www.anton-paar.com/ca-en/products/details/rheome>
Cup and Bob:
<http://www.anton-paar.com/ca-en/products/details/rheome>
<https://ciks.cbt.nist.gov/~garbocz/SP946/node14.htm>

SANS Measurement




Sample Preparation

NaCl Concentrations

30 mM

150 mM

450 mM



0.6 mM

1.2 mM

0.6 mM

1.2 mM

1.8 mM

0.6 mM


1.2 mM

1.8 mM



Insulin Concentrations

- Insulin
- NaCl
- Zinc
- D₂O
- Phenol
- Phosphate Buffer


NaCl ions screening charges



Rheology: Viscosity vs. Shear Stress




Insulin Concentration	Volume Percentage
0.6mM	0.44%
1.2mM	0.88%
1.8mM	1.32%




SANS Data: Increasing Salt Concentration


0.6mM Insulin




Radial Configuration



Tangential Configuration





Isotropic Orientation




Higher Order Structures: SAXS Comparison

450mM NaCl




Shear Thickening Dependence on Shape


0.6mM Insulin




30mM NaCl



150mM NaCl






450mM NaCl



SANS Data: Increasing Salt Concentration


--- Fitted models




**Shear-induced reversible structural change at 150mM NaCl, 3500s⁻¹

Electrostatic Interactions?



30mM NaCl



150mM NaCl




450mM NaCl




SHEAR THICKENING

Electrostatic Interactions?


30mM NaCl




150mM NaCl




450mM NaCl



With NaCl Screening



Easier to overcome repulsion



Higher NaCl = Lower σ^*

Summary of Findings

Structure

Rods **form, lengthen, and order** with increasing NaCl and insulin concentrations

Very little shear dependence

Rheology

Shear-thinning \longrightarrow Shear-thickening

Critical shear stress (electrostatic barrier?)

Future Outlook

Tests with higher shear rates, lower q

More rheological studies

Acknowledgments

Center for High Resolution Neutron Scattering

Novo Nordisk

Grethe Vestergaard Jensen, PhD

Mentor, NIST Center for Neutron Research

Joe Dura, PhD and Julie Borchers, PhD

NCNR SURF Coordinators, NIST Center for Neutron Research

Brandi Toliver, PhD

SURF Program Director, National Institute of Standards and Technology

